

April 4th, 2024
Security Notice CVE-2024-3094

XZ LZMA Vulnerability

Dear Sir or Madam,

KARL STORZ is aware of and currently monitoring the XZ LZMA Vulnerability (CVE-2024-3094). This vulnerability was announced on March 29th, 2024. For the moment we can confirm that KARL STORZ or our customers are not harmed in any way. Certainly, we will keep monitoring and investigating that important matter and will keep you informed.

For the moment, we can assure you that all KARL STORZ medical products, non-medical product or their accompanying infrastructure do not utilize the affected version 5.6.0 and 5.6.1 of the library liblzma.

KARL STORZ always recommends the following steps:

Employ Good Network Hygiene Practices

- Ensure data has been backed up and stored according to your individual processes and disasterrecovery procedures.
- Execute updates to malware protection, where available.

For more information on this CVE, please follow the link below:

<https://nvd.nist.gov/vuln/detail/CVE-2024-3094>

If you have questions, please contact your KARL STORZ sales representative, or contact us at the following address: productsecurity@karlstorz.com.

Revision History

Version	Date	Comment
1.0	2024-03-30	Initial advisory

Annex I - List of products

Article/Model	Product Description	Software Version	Affected by CVE-2024-3094
All	All KARL STORZ products	All software versions	Not Affected